

Variolink[®] Veneer

ESTHETIC RESIN CEMENT

ivoclar
vivadent[®]

EXCEPTIONAL ESTHETICS

Value Shading

Variolink Veneer offers seven light cure value shades to enhance the esthetics of the final restoration. The Medium Value (MV 0) is very translucent and has little effect on the final restoration. The High Value shades allow the clinician to brighten or increase the value and the Low Value shades provide a "warming" effect to the final restoration.



Variolink Veneer Light Cure

All Variolink Veneer shades are amine-reduced providing a high degree of color stability ensuring outstanding esthetics.



Variolink[®] Veneer

EXTRAORDINARY CLEAN UP

Clean Up

Variolink Veneer offers the clinician a fast, easy solution to cleaning up excess cement during the adhesive bonding procedure. No more grinding, burs, or tedious clean up required. Simply load the restorations, seat, and tack & wave with a dental curing light. Excess cement can be easily peeled away. Prior to complete polymerization Variolink Veneer provides the unique combination of the security of adhesive bonding and the easy clean up of a conventional cement.




Clinical photography courtesy of David S. Hornbrook, DDS, FAACD
Executive Director, Hornbrook Group

Physical Properties


Variolink Veneer's microfill chemistry provides key benefits including low wear rates, high translucency, high polishability and strong bonds. Variolink Veneer demonstrates significantly lower wear values compared to competitive resin cements currently on the market ensuring the long term clinical success of your esthetic restorations.

TRANSLUCENCY



* Not registered trademarks of Ivoclar Vivadent, Inc.
2003, R&D Ivoclar Vivadent, Inc.

WEAR



* No registered trademarks of Ivoclar Vivadent, Inc.
UAB School of Dentistry, October 2003
Shiro Suzuki, DDS, PhD, FAD

Variolink® Veneer

DELIVERY FORMS



597614 Variolink Veneer Assortment

- Contents:**
- 7 syringes Variolink Veneer, 1.5 ml each, in 7 shades
 - 7 syringes Variolink Veneer Try-In, 1.5 ml each, in 7 shades
 - 1 bottle Monobond S, 5 g
 - 1 syringe Liquid Strip, 2.5 g
 - 1 syringe Total Etch, 2 g
 - 12 Excite single dose vessels, 0.1 g each
 - 1 syringe Total Etch, 2 g
 - Various accessories



Variolink Veneer Refills

Contents: 1 syringe Variolink Veneer, 1.5 ml, available in the following shades:

591423	High Value	+3
591424	High Value	+2
591425	High Value	+1
591426	Medium Value	0
591427	Low Value	-1
591428	Low Value	-2
591429	Low Value	-3



Variolink Veneer Try-In Refills

Contents: 1 syringe Variolink Veneer Try-In, 1.5 ml, available in the following shades:

591430	High Value	+3
591431	High Value	+2
591432	High Value	+1
591433	Medium Value	0
591434	Low Value	-1
591435	Low Value	-2
591436	Low Value	-3

100% CUSTOMER SATISFACTION
GUARANTEED!

UNITED STATES

175 Pineview Drive
Amherst, NY 14228
716-691-0010
800-533-6825
Fax: 716-691-2285
www.ivoclarvivadent.us.

CANADA

2785 Skymark Ave., Unit 1
Mississauga, ON L4W 4Y3
(905) 238-5700
(800) 263-8182
Fax: 905-238-5711
www.ivoclarvivadent.ca

